

# Hematocrit Interference is a Common Phenomenon in Many Devices for Glucose Self-Measurement

Ramljak S<sup>1</sup>, Lock JP<sup>2</sup>, Schipper C<sup>1</sup>, Musholt PB<sup>1</sup>, Forst T<sup>1</sup>, Pfützner A<sup>1</sup>

<sup>1</sup>IKFE – Institut für klinische Forschung und Entwicklung GmbH, Mainz  
<sup>2</sup>University of Massachusetts, Worcester, MA

## Background


High hematocrit levels may lead to false low glucose readings and vice versa. Aim of this laboratory investigation was to assess the potential influence of hematocrit variations on a variety of blood glucose meters applying different measurement technologies.

## Methods



Venous heparinized blood was manipulated to contain 3 different blood glucose concentrations (50-90 mg/dL, 120-180 mg/dL, and 280-350 mg/dL) and 5 different HCT values (25%, 35%, 45%, 55%, and 65%). After careful oxygenation to capillary blood oxygen pressure, each sample was measured eight times with the following devices: NovaMax Link & NovaMax+, AccuChek Aviva, Nano & Active, BGStar, Contour & Breeze, OneTouch Ultra2 & Verio2, Freestyle Freedom Lite, Precision Xtra & Optium Xceed, GlucoMen LX+ GlucoFix Mio+, OnCall+ & Platinum, GlucoLab, GlucoDr, and TaiDoc fora. YSI 2300 served as plasma reference method. Stability to hematocrit influence was assumed, with <10 % mean glucose result deviation between the highest and lowest hematocrit levels.

## Results

Seven of the investigated meters showed a stable performance in this investigation: BG\*Star (3%), Verio2 (3%), NovaMax+ (4%), GlucoFix Mio+ (4%), GlucoMen LX+ (5%), NovaMax Link (7%), AccuChek active (7%). All other meters failed this hematocrit interference test, with FreeStyle Freedom lite (11%), Platinum (12%) being the better devices and OnCall+ (68%), GlucoLab (51%), TaiDoc (39%), Breeze (38%) showing the worst performance.



HIF describes the sum of the maximal deviation of the mean HCT curve above 100 % plus the maximal deviation of the mean HCT curve below the 100 % reference (see graph).



Based on clinical and laboratory considerations, a HIF factor below 10 % has been arbitrarily defined as acceptable

## Conclusion

Hematocrit variations may occur frequently in daily routine (e.g., due to dehydration/exercise, nicotin and alcohol abuse, pregnancy etc.). Our results encourage use of meters with stable performance under these conditions.